


版本 REV.	描述 DESCRIPTION	日期 DATE	批准 APPD
A	新制成 New made	2023.03.29	GJM



备注：此配件单独包装,置于内包装袋中。

P8	P1XS	止液夹 Clamp, WALL1/8"-WALL1/4"	4 pcs
P7	P32	扎带 Cable Tie, 2.5 mm	4 pcs
P6	PCG25	支撑杆 Support Rod, 250 mm	2 pcs
P5	P41+P51	防尘袋 Dust Cover	2 pcs
P4	P31	扎带 Cable Tie, 3.6 mm	A/R
P3	P3T185	StepLess® 卡箍 Ear Clamps, 18.5 mm (仅供组装参考 For assembly reference only)	18 pcs
P2	P1RH	止液夹 Clamp, OD1/4"-OD1/2"	4 pcs
P1	P1HM	止液夹 Clamp, OD1/2"-OD3/4"	6 pcs
7	P8S	快装盲盖 TC Blind Cap, 1-1/2"TC	3 pcs
6	P9S	快装卡箍 TC Clamp, 1-1/2"TC	5 pcs
5	PBS	硅胶密封圈 Silicone Gasket, 1-1/2"TC	2 pcs
4	C4YS	宝塔转快装接头 TC Fitting, 3/8"HB to 1-1/2"TC	2 pcs
3	SFGSY	仪表三通接头 Instrument Tee, 3/8"HB-3/8"HB-1/2"TC, 焊膜	3 pcs
2	C6TYY	T型等径三通接头 T-Fitting, 3/8"HB	5 pcs
1	AQGI7006	AseptiQuik® G 无菌连接器 Connector, 3/8"HB	2 pcs
F4	C1FYV1	Bricap® C01 Capsule Filter with Cobetter Corevital® SMD 0.2 μm, 货号 Cat. No.: C01FYSDNG1P	1 pcs
F3	V5DBB ACEBG1P	Bricap® V50D Capsule Filter with Cobetter Aegivast® 0.2 μm, 货号 Cat. No.: V5DBBACEBG1P	2 pcs
F2	C2TYA2	Bricap® C02 Capsule Filter with Cobetter Corevital® APB 0.2 μm, 货号 Cat. No.: C02TYAPBBG1P	2 pcs
F1	C3FYE1	Bricap® C03 Capsule Filter with Cobetter Purcise® SLE 0.2 μm, 货号 Cat. No.: C03FYSLENG1P	2 pcs
T3	TVY	Vena® Technosil 编织增强铂金硫化硅胶管 Reinforced Silicone Tubing, ID9.6*OD17.0 mm	110 cm
T2	TP15	Pharma-50 铂金硫化硅胶管 Pt-cured Silicone Tubing, ID3/16*OD3/8"	60 cm
T1	TPY	Pharma-50 铂金硫化硅胶管 Platinum Cured Silicone Tubing, ID3/8*OD5/8"	230 cm
A2	K-2D050 CBRR	500 ml 1/4"-1/4" 接口储液袋 500 mL Single-use Bag with 1/4"-1/4" Ports	2 pcs
A1	K-2D05L C3YNY	5 L 3/8"-3/8" 接口储液袋 5 L Single-use Bag with 3/8"-3/8" Ports	2 pcs
序号 NO.	型号 PART NO.	名称 DESCRIPTION	数量 QTY

γ 灭菌 IRRADIATED 25-45 kGy	材料 MATERIAL NA	 杭州科百特过滤器材有限公司 www.cobetterfiltration.com		
制图 DRN ZHX	审核 CKD ZJH			批准 APPD GJM
2023.03.29	2023.03.29			2023.03.29
比例 SCALE NTS	版本 REV. A	页码 SHEET 1 of 1	产品货号 CAT. NO. SAT6XY40AC3E1AS01	
图纸编号 DRAW'G NO. C110-23070	图纸名称 TITLE Lifecube® SA 一次性使用囊式过滤器组件 Single-use Assemblies with Capsule Filter			

必填 * 请在应用点勾选 "√"
Required* Please check '√' in the application.

上游产品 Upstream 下游产品 Downstream 原液 Drug Substance 制剂 Drug Product 其他 Others _____

客户确认：Customer Approval 日期：Date

该图纸包含杭州科百特过滤器材有限公司专有的保密信息，未经本公司事先书面授权，不得复制、公开或以其他方式使用该图纸。
The confidential information contained within this drawing is proprietary to Cobetter, and shall not be reproduced, disclosed, or used in other ways without prior written authorize of Cobetter Filtration Equipment Co., Ltd.

© 2023 Cobetter
保密文件 CONFIDENTIAL