


版本 REV.	描述 DESCRIPTION	日期 DATE	批准 APPD
A	新制成 New made	2023.02.09	GJM



P2	P31	扎带 Cable Tie, 3.6 mm	A/R
P1	P1RH	止液夹 Clamp, OD1/4"-OD1/2"	27 pcs
5	C52	鲁尔螺帽 Luer Thread Cap	10 pcs
4	C51	无针取样器 Needlefree Swabable Valve	10 pcs
3	C4RC	宝塔转母鲁尔接头 1/4"HB to Female Luer Lock	10 pcs
2	C6XRRRR	四通转接头 Cross Fitting, 1/4"HB	5 pcs
1	PAR	软管堵头 Tubing Plug, 1/4"HB	1 pcs
F1	V5DBB ACEBG1P	Bricap® V50D Capsule Filter with Cobetter Aegivast® 0.2 μm, 货号 Cat. No.: V5DBBACEBG1P	1 pcs
T1	TCR	C-Flex® 374 热塑管 TPE Tubing, ID1/4"OD7/16"	420 cm
A2	K-2D010 CBRR	100 ml 1/4"-1/4" 接口储液袋 100 mL Single-use Bag with 1/4"-1/4" Ports	10 pcs
A1	K-2D02L C3RNR	2 L 1/4"-1/4" 接口储液袋 2 L Single-use Bag with 1/4"-1/4" Ports	1 pcs
序号 NO.	型号 PART NO.	名称 DESCRIPTION	数量 QTY

γ 灭菌 IRRADIATED 25-45 kGy		材料 MATERIAL NA		 杭州科百特过滤器材有限公司 www.cobetterfiltration.com
制图 DRN ZHX 2023.02.09	审核 CKD ZJH 2023.02.09	批准 APPD GJM 2023.02.09		
比例 SCALE NTS	版本 REV. A	页码 SHEET 1 of 1	产品货号 CAT. NO. KAS010XCPF01	
图纸编号 DRAW'G NO. C110-23049		图纸名称 TITLE Lifecube® AS 一次性使用取样组件 Single-use Sampling Assemblies		

必填 \* 请在应用点勾选 "√"  
Required\* Please check '√'  
in the application.

上游产品 Upstream  下游产品 Downstream  原液 Drug Substance  制剂 Drug Product  其他 Others \_\_\_\_\_

客户确认：  
Customer Approval

日期：  
Date

该图纸包含杭州科百特过滤器材有限公司专有的保密信息，未经本公司事先书面授权，不得复制、公开或以其他方式使用该图纸。  
The confidential information contained within this drawing is proprietary to Cobetter, and shall not be reproduced, disclosed, or used in other ways without prior written authorize of Cobetter Filtration Equipment Co., Ltd.

© 2023 Cobetter  
保密文件 CONFIDENTIAL